

AFM Data Processing made easy: The latest Bruker BioAFM developments


Alexander Dulebo

Bruker Nano GmbH, JPK BioAFM, Am Studio 2D, 12489 Berlin, Germany


We will demonstrate some new and useful features of the AFM data processing software, developed by the Bruker BioAFM software team: **NanoScope Analysis** and **Data Processing**. The intuitive software enables users generate high-quality, correlative data.

Among the topics covered, we will briefly discuss:

- Advance cropping of force curve maps (QI)
- Batch data analysis of multiple folders
- Tracing analysis process for every file
- Nano Mechanics: JKR and Cone Sphere models
- Creating AFM-optical overlays offline
- And more



Height (measured) (qi-fit-2019.09.05)



Microscope Image Registration and Overlay (MIRO)

MIRO offline enables import, registration and overlays of AFM images and datasets, additionally enabling selection of datasets for further analysis.

